


Product Data Sheet

AQUA KEM® BLUE LAVENDER CONCENTRATED



Product	Aqua Kem® Blue Lavender Concentrated													
Short product description	Concentrated toilet fluid for the waste-holding tank													
Consumer proposition	A powerful, highly effective, 2½ times concentrated toilet fluid with a fresh lavender fragrance for the waste-holding tank of your Cassette and portable toilet.													
Volume (depends per country)	0.78 L													
Color	Blue													
Perfume	Lavender													
Way to use	<p>Dosage for a waste-holding tank of up to 20 L:</p> <table border="1"> <thead> <tr> <th>Ambient temperature</th> <th>Dose</th> <th>Fresh water</th> </tr> </thead> <tbody> <tr> <td>> 30°C</td> <td>70 ml</td> <td>2-3 l</td> </tr> <tr> <td>Standard dose</td> <td>60 ml</td> <td>2-3 l</td> </tr> <tr> <td>< 15°C</td> <td>50 ml</td> <td>2-3 l</td> </tr> </tbody> </table> <p>1 dosage lasts up to 5 days.</p>		Ambient temperature	Dose	Fresh water	> 30°C	70 ml	2-3 l	Standard dose	60 ml	2-3 l	< 15°C	50 ml	2-3 l
Ambient temperature	Dose	Fresh water												
> 30°C	70 ml	2-3 l												
Standard dose	60 ml	2-3 l												
< 15°C	50 ml	2-3 l												
General USP's	<ol style="list-style-type: none"> 0.78L Aqua Kem® Blue Lavender CONCENTRATED = 2L Aqua Kem® Blue Lavender Smaller bottle, bigger impact* Significantly improved waste liquefaction for easier tank emptying* Significantly improved reduction of gas build-up* Equal outstanding odour control* Equal fresh lavender fragrance* Equal number of dosages* Less weight* Easy to handle bottle 													

	*) compared to the 2L Aqua Kem® Blue Lavender
Positioning	In Thetford's concentrated Aqua range of additives
Performance	1 dosage lasts up to 5 days
EU/GHS - CLP classification	 <p>Danger Contains: 2-BROMO-2-NITROPROPANE-1,3-DIOL (H318) Causes serious eye damage (H400) Very toxic to aquatic life (P101) If medical advice is needed, have product container or label at hand (P102) Keep out of reach of children (P280) Wear eye protection (P305+P351+P338+P310) IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician (P391) Collect spillage (P501) Dispose of contents/container to a hazardous or special waste collection point</p>
Transport classification (ADR)	UN-no. 3082 (LQ)
Shelf life	3 years. After the 'best used by date', which is printed on the bottle, the additive may still be used. There may be a slight reduction in performance.
SDS	Available through wholesaler
Available pictures	Pictures available on website (Low and High resolution)

Order information				
Available in countries	UK / Germany / France / Benelux / Italy / Spain / Portugal / Denmark / Norway / Sweden / Finland / Poland / Greece / Russia / China / Australia			
Weight bottle + liquid	878,5 g			
Quantity outer case	12 x 0.78L			
Measurements outer case	392 x 194 x 233			
Quantity on pallet	72 boxes			
	Languages	Countries	Part number box	EAN-code bottle
Language combination 1	English / German / Polish	UK/GER/POL	30625CX	8710315025996
Language combination 2	Dutch / French / German	BENLX/FRA	30625AA	8710315025965
Language combination 3	Italian / Spanish / Portuguese	ITA/SPA/PORT	30625CQ	8710315025958
Language combination 4	Danish / Norwegian / Russian	DEN/NOR/RUS	30625CW	8710315025989
Language combination 5	Swedish / Finnish / Greek	SWE/FIN/GRE	30625CV	8710315025972
Language combination 6	English / Chinese	CHINA only	30625CK	8710315026009
Language combination 7	English	AUSTALIA only	30625ZK	8710315026016