



Product Data Sheet

AQUA KEM® GREEN CONCENTRATED



Product	Aqua Kem® Green Concentrated												
Short product description	Concentrated toilet fluid for the waste-holding tank												
Consumer proposition	Environmentally conscious, highly effective, 2 times concentrated toilet fluid for the waste-holding tank of your Cassette and portable toilet.												
Volume (depends per country)	0.75 L												
Color	Green												
Perfume	Citrus fruits												
Way to use	<p>Dosage for a waste-holding tank of up to 20 L:</p> <table border="1"> <thead> <tr> <th>Ambient temperature</th> <th>Dose</th> <th>Fresh water</th> </tr> </thead> <tbody> <tr> <td>> 30 °C</td> <td>85 ml</td> <td>2-3 L</td> </tr> <tr> <td>Standard dose</td> <td>75 ml</td> <td>2-3 L</td> </tr> <tr> <td>< 15 °C</td> <td>65 ml</td> <td>2-3 L</td> </tr> </tbody> </table> <p><i>1 dosage lasts up to 4 days.</i></p>	Ambient temperature	Dose	Fresh water	> 30 °C	85 ml	2-3 L	Standard dose	75 ml	2-3 L	< 15 °C	65 ml	2-3 L
Ambient temperature	Dose	Fresh water											
> 30 °C	85 ml	2-3 L											
Standard dose	75 ml	2-3 L											
< 15 °C	65 ml	2-3 L											
General USP's	<ol style="list-style-type: none"> 0.75L Aqua Kem® Green CONCENTRATED = 1.5L Aqua Kem® Green Smaller bottle, bigger impact* 												

	<ul style="list-style-type: none"> 3. Septic Tank Safe (Thetford anaerobic biodegradability test) 4. Significantly improved waste liquefaction for easier tank emptying* 5. Significantly improved reduction of gas build-up* 6. Equal outstanding odour control* 7. Equal fragrance* 8. Equal number of dosages* 9. Less weight* 10. Easy to handle bottle 11. Keeps the waste-holding tank even cleaner than before <p><i>*) compared to The Original 1.5L Aqua Kem® Green</i></p>
Positioning	In Thetford's concentrated Aqua range of additives
Performance	1 dosage lasts up to 4 days
EU/GHS - CLP classification	 <p>Danger Contains: 2-BROMO-2-NITROPROPANE-1,3-DIOL; calcium dinitrate (H318) Causes serious eye damage. (P101) If medical advice is needed, have product container or label at hand. (P102) Keep out of reach of children. (P103) Read label before use. (P305+P351+P338+P310) IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a doctor.</p>
Transport classification (ADR)	none
Shelf life	3 years. After the 'best used by date', which is printed on the bottle, the additive may still be used. There may be a slight reduction in performance.
SDS	Available through wholesaler
Available pictures	Pictures available on website (Low and High resolution)

Order information				
Available in countries	UK / Germany / France / Benelux / Italy / Spain / Portugal / Denmark / Norway / Sweden / Finland / Poland / Russia / Australia			
Weight bottle + liquid	916 g			
Quantity outer case	12 x 0.75L			
Measurements outer case	392 x 194 x 233			
Quantity on pallet	72 boxes			
	Languages	Countries	Part number box	EAN-code bottle
Language combination 1	English / German / Polish	UK/GER/POL	30646CX	8710315995220
Language combination 2	Dutch / French / German	BENLX/FRA	30646AA	8710315995237
Language combination 3	Italian / Spanish / Portuguese	ITA/SPA/PORT	30646CQ	8710315995244
Language combination 4	Danish / Norwegian / Russian	DEN/NOR/RUS	30646CW	8710315995251
Language combination 5	Swedish / Finnish	SWE/FIN	30646AF	8710315995268
Language combination 7	English	AUSTALIA only	30646ZK	8710315995275